

Foam Tipped Applicator

Polypropylene Shaft

1876-PF / 1876-PFB

NVR

NON-VOLATILE RESIDUE

Test method IDEMA method M7-98

NVR Solvent - Hexane

NVR 488 µg/swab

The minimum detection limit for NVR is 10 µg/swab

ION CHROMATOGRAPHY TESTING RESULTS

IONIC CLEANLINESS

Test method IDEMA M13-99 (anions) M12-99 (cations)

ANIONS

Fluoride	ND	µg/swab
Chloride	0.6320	µg/swab
Phosphate	ND	µg/swab

CATIONS

Lithium	ND	µg/swab
Sodium	0.1080	µg/swab
Ammonium	0.0280	µg/swab
Potassium	0.0600	µg/swab
Magnesium	ND	µg/swab
Calcium	0.0960	µg/swab

ND - Non Detected

Test procedures available upon request.



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

To Order Call:
Harmony Business Supplies
1-800-899-1255
www.HarmonyCR.com