

1616-PF PurSwab Foam Tipped Cleaning Applicator



Features and Benefits

Tip: • Open cell 100 PPI polyurethane foam • 6.35mm diameter • 22.2mm length	Shaft: • Extruded white polypropylene • 3mm diameter • 150.8mm overall length	Bond: • Thermal	Uses: • Cleaning • Applying adhesives	Packaging: • 20 bags of 50 each • Resealable poly bags	Qualities: • Soft, lint and residue free • Can be used with IPA or Deionized water, or dry • Thickwalled sturdy shaft
---	--	--	---	--	--



Specialty Applicators

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

To Order Call:
Harmony Business Supplies
1-800-899-1255
www.HarmonyCR.com