

# 1803-PF PurSwab Foam Tipped Applicator



## Features and Benefits

<b>Tip:</b> <ul style="list-style-type: none"> <li>• Reticulated 100PPI medical grade foam</li> <li>• 11.88mm length</li> <li>• 3.96mm width</li> <li>• Mini flexible paddle tip</li> <li>• Bonded securely to shaft with thermo-plastic adhesive</li> </ul>	<b>Shaft:</b> <ul style="list-style-type: none"> <li>• Black polypropylene</li> <li>• 76.14mm length</li> <li>• 2.31mm diameter</li> </ul>	<b>Standards:</b> <ul style="list-style-type: none"> <li>• ISO 9001:2000</li> <li>• Produced under the FDA's GMP's</li> </ul>	<b>Uses:</b> <ul style="list-style-type: none"> <li>• Typically used with deionized water and isopropyl</li> <li>• Removing epoxy residue from connectors</li> <li>• General cleaning</li> <li>• Flux remover</li> </ul>	<b>Packaging:</b> <ul style="list-style-type: none"> <li>• 1,000 per case (20 bags of 50 in resealable poly zip bags or 2 bags of 500)</li> </ul>	<b>Qualities:</b> <ul style="list-style-type: none"> <li>• Soft non-abrasive foam</li> <li>• Can be used with IPA and DI water, or dry</li> <li>• Lint and residue free</li> </ul>
--	--	---	--	---	--



Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

**Tel** 800-321-2313 (US and Canada Only)  
207-876-3311

**Fax** 800-323-4153  
207-876-3130

**Web** [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

**Email** [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

To Order Call:  
**Harmony Business Supplies**  
1-800-899-1255  
[www.HarmonyCR.com](http://www.HarmonyCR.com)