

3677 PurSwab Microfiber Tipped Applicator



Features and Benefits

Tip:	Shaft:	Uses:	Packaging:	Qualities:	Extraction Analysis*
• Ultra clean microfiber • Flexible paddle • 6.8mm width • 16.8mm length • Ultra sonic welded	• Molded polypropylene • Solid • 151.7mm length • 3mm diameter	• Cleaning electronics and equipment • Removing residues and contaminants • Applying solution	• 10 bags of 100 • Double poly bags	• Superior wipe dry characteristics • Low in particulates • Soft non-abrasive • NVR 0.003 mg/swab • Lint free • Can be used dry or with IPA, di water or solvents	• Chloride (Cl) 0.28 µg/tip • Nitrite 0.275 µg/tip • Bromide <0.025 µg/tip • Nitrate 0.15 µg/tip • Phosphate 0.515 µg/tip • Sulphate 0.185 µg/tip • Sodium 0.312 µg/tip • Ammonium 0.03 µg/tip • Potassium 0.051 µg/tip • Magnesium 0.723 µg/tip • Calcium 2.955 µg/tip
					* Based on a typical sample batch



Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

To Order Call:
Harmony Business Supplies
1-800-899-1255
www.HarmonyCR.com