

HEEL GRINDER ESTAHEEL

The heel grinder Estaheel is used on the heel of the shoe wear, which through a conductive strap and a rubber part, puts the human body in contact with the ground to allow electrostatic charges to go through.

> Benefits


- Traceability through bar code
- Longer lifetime due to its encapsulated resistance
- Excellent abrasion resistance
- Does not damage shoe wear
- Made in Mexico



EHL-1001

EXCELLENT DURABILITY RESISTANCE

> General information

Size / Length	Unitalla
UM	Pieces
Production / Sale batch	10 Pieces
Production / Sale batch	400 Pieces
Colors	Royal blue and yellow 

Industry	Electronics, automotive and related
Cleaning	Use a damp cloth

> Technical Information

Lower limit (ohms)	$8 \times 10 e^{-5}$
Upper limit (ohms)	$3.5 \times 10 e^{-7}$
Voltage generation	N / A
Download time	0.01 sec

> Standards , warranties and care

Rules applying	ANSI/ESD S20.20-2014 ANSI/ESD SP9.2-2003
Physical guarantee	6 months
Electric guarantee	6 months
*Physical guarantee	6 months
*Electric guarantee	3 months

*With unique heel grinder

> Contradictions

Do not expose to chemical agents, aggressive substances, and highly abrasive surfaces such as pavement, stone lined paths, dirt roads, etc.

These data are based on information we believe correct; are offered in good faith but without guarantee in respect of the application and use of the product, since the conditions and methods of use are beyond control of the company. Before the final use of the product, we recommend the user to perform a detailed assessment of the same through the implementation and use of significant samples. Users are responsible for checking the validity of the information provided.

ESTATEC UNIQUE TRACKING CODE (CURE)

The Estatec unique tracking code is an identification code that helps on product tracking. It's used to know data from production dates as well as date of delivery and guarantees. It also has a bar code for easy reading using electronic methods.

In the heel grounder, a part of the label stays encapsulated so that there is always a backup in case of damage to the exterior tag or label



CURE EXAMPLE

MODELS CHART

> Code	> Name	> Description	> Production batch pieces	> Standard packing pieces
EHL-1001	Estaheel talonera azul 1 MΩ	standard model, 2.5 cm wide. with conductive tape	10	300
EHL-1002	Estaheel talonera azul 2 MΩ	Resistance of 2 milliohms	100	300
EHL-1041	Estaheel talonera amarilla 1 M	Classic model , Yellow , Resistor 1 Mohm , Lamina Width 2.5 Cm , Unitalla	100	300
EHL-1061	Estaheel talonera amarilla 1 MΩ	One sole side , 2.5 cm wide, elastic back	100	300
EHL-1007	Talonera con azul rey 1 MΩ	Model long conductive tape (60cm) , royal blue , resistance 1 M , Lamina 2.5 width , size fits all	100	300

These data are based on information we believe correct; are offered in good faith but without guarantee in respect of the application and use of the product, since the conditions and methods of use are beyond control of the company. Before the final use of the product, we recommend the user to perform a detailed assessment of the same through the implementation and use of significant samples. Users are responsible for checking the validity of the information provided.