

SOLE GRINDER

The sole , It is used in the heel footwear, which , by tape and rubber, puts the human body in contact with the ground to allow the passage of electrostatic charges.

> Benefits


- Traceability through bar code
- Longer lifetime due to its encapsulated resistance
- Excellent abrasion resistance
- Does not damage shoe wear
- Made in Mexico



EHL-1023

HAS BIGGER CONTACT AREA AND GREATER TEAR RESISTANCE

> General information

Size / Length	S, M & L
UM	Pieces
Production / Sale batch	10 Pieces
Colors	Royal blue 

Industry	Electronics, automotive and related
Cleaning	Use a damp cloth

> Technical Information

Lower limit (ohms)	2×10^6
Upper limit (ohms)	3.5×10^7
Voltage generation	N / A
Download time	0.01 sec

> Standards , warranties and care

Rules applying	ANSI/ESD S20.20-2014 ESD SP9.2-2003 TR53-01-06
Physical guarantee	6 months
Electric guarantee	6 months

> Contradictions

Do not expose to chemical agents, aggressive substances, and highly abrasive surfaces such as pavement, stone lined paths, dirt roads, etc.

These data are based on information we believe correct; are offered in good faith but without guarantee in respect of the application and use of the product, since the conditions and methods of use are beyond control of the company. Before the final use of the product, we recommend the user to perform a detailed assessment of the same through the implementation and use of significant samples. Users are responsible for checking the validity of the information provided.

ESTATEC UNIQUE TRACKING CODE (CURE)

The Estatec unique tracking code is an identification code that helps on product tracking. It's used to know data from production dates as well as date of delivery and guarantees. It also has a bar code for easy reading using electronic methods.

In the sole case, a part of the label stays encapsulated so that there is always a backup in case of damage to the exterior tag or label



CURE EXAMPLE

MODELS CHART

> Code	> Name	> Description	> Production batch pieces	> Standard packing pieces
EHL-1022	Blue rubber sole, 2 MΩ -S-	2 MΩ resistance reinforced rubber on sole and heel. Size S	10	300
EHL-1023	Blue rubber sole, 2 MΩ -M-	2 MΩ resistance reinforced rubber On sole and heel. Size M	10	300
EHL-1024	Blue rubber sole, 2 MΩ -L-	2 MΩ resistance reinforced rubber On sole and heel. Size L	10	300

These data are based on information we believe correct; are offered in good faith but without guarantee in respect of the application and use of the product, since the conditions and methods of use are beyond control of the company. Before the final use of the product, we recommend the user to perform a detailed assessment of the same through the implementation and use of significant samples. Users are responsible for checking the validity of the information provided.